

Valérie Lemieux

Vice Présidente Ventes et Opérations

514-337-1184 / Fax: 514-337-0913

valerie@dardo.ca

www.dardo.ca



DARDO
ORTHOPÉDIQUE
ORTHOPÉDIE



extrostyle
medical shoes



CERTIFICATIONS & RELIABILITY

All ExtroStyle footwear has the CE mark and is certified as Class I Medical Devices (in compliance with MDR 745/17), registered at the Italian Ministry of Health. Our products also meet the requirements of UNI EN ISO 9001 and UNI EN ISO 13485 standards.



RECOMMENDED PRODUCT BY RISK PHASE

LEVEL OF DIABETIC RISK	PROBLEM	RECOMMENDED PRODUCT
0 LOW RISK	Diabetic patient with normal conformed foot. No neuropathy	Energy Line Comfort Line Rehabilitation Line
1 MEDIUM RISK	Neuropathy with light deformities	Energy Line Rehabilitation Line
2 HIGH RISK	Neuropathy or arterial disease. Relevant deformities but without ulceration and/or previous amputation.	Energy Line Rehabilitation Line
3 VERY HIGH RISK	Neuropathy or arterial disease. Deformities, ulceration and/or previous amputation.	Energy Line Rehabilitation Line

ENERGY LINE

The reinforced, solid heel support surrounds and protects the heel, preventing rubbing and injuries.

The padded, comfortable, and breathable upper—designed without seams—envelops the foot in a protective and cushioning shell.

Designed to accommodate custom mainsoles up to 9 mm thickness


Breathable upper

Roller-profile outsole with front and rear curvature

Improved weight distribution that promotes a smooth and safe gait

The insoles used for the ENERGY lines are composed of a high-density open-cell foam that ensures absorption and deabsorption capacity, keeping the foot ventilated, dry and odourless, as well as ensuring a pleasant sensation when walking.

The ergonomic design of the footbed offers high comfort, minimising the effort and tensions in the plantar support, both during walking and in a static position, dampening the impact of pressure during walking.



CHARACTERISTICS OF THE ENERGY LINE INSOLES

Insoles are available in 3 different thickness:

- ① 9 mm footbed (by default)
- ② 6 mm footbed
- ③ 3 mm footbed



①

②

③

LINE ENERGY
EXTRA WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



SUITABLE FOR THE
PREVENTION OF
DIABETIC FOOT



PADDED
LINING

ENERGY MARINA
Certification no. 2702692

Diabetic foot risk level
0 1 2 3 See page 2



Color Black
Sizes 35 to 43



ADJUSTABLE DOUBLE
VELCRO OPENING

ENERGY LINDA
Certification no. 2702663

Diabetic foot risk level
0 1 2 3 See page 2



ADJUSTABLE DOUBLE
VELCRO OPENING

ENERGY ELENA
Certification no. 2702681

Diabetic foot risk level
0 1 2 3 See page 2



Color Black
Sizes 35 to 43



LINE ENERGY
EXTRA WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



SUITABLE FOR THE
PREVENTION OF
DIABETIC FOOT



PADDED
LINING

ENERGY LUIGI

Certification no. 2702694

Diabetic foot risk level

0 1 2 3

See page 2



Color Black

Sizes 39 to 50



BIOMECHANICAL
OUTSOLE

ENERGY ALEX NERO

Certification no. 2702691

Diabetic foot risk level

0 1 2 3

See page 2



Color Black

Sizes 39 to 50

ENERGY DANIELE

Certification no. 2702672

Diabetic foot risk level

0 1 2 3

See page 2



Color Black

Sizes 35 to 50




TOTAL OPENING WITH
DOUBLE ZIP

PHASE1 PREVENTION

Thanks to the predisposition for the Energy 10 **removable insole**, it is a great ally for the prevention or containment of the most common foot diseases such as: hallux valgus, hammertoes, metatarsalgia, heelitis and vascular diseases as well as for the prevention with a 0 - 1 risk in diabetic foot.

The **flexible biomechanical energy outsole** allows a natural transition between the three phases of

walking, ensuring a reduction of the pressure peaks on the forefoot and limiting the wear on the joints, thus allowing a more spontaneous walk and reducing the feeling of fatigue. The **self-moulding and seamless upper** at the rubbing points envelops the anatomy of the foot and its possible deformities, ensuring comfort and preventing painful friction.



DIABETIC FOOT
VASCULOPATHIES



HALLUX VALGUS



HAMMER TOES
CLAW TOES
METATARSALGIA



RHEUMATIC FOOT
MORTON'S NEUROMA



BURSITIS
TARSAL TUNNEL



CHARCOT FOOT

PHASE2 POST-OPERATIVE REHABILITATION

The energy 10 full opening shoe with 2 zips, equipped with a flexible outsole, can be defined as a multifunctional shoe thanks to a combination of insoles that are variable in height. **It can be worn from the first day after the operation until the end of convalescence.**

This implies a more correct and balanced walking and distribution of body weight, as well as an undeniable benefit for the back, pelvis and legs. Also the psychological well-being is not to be underestimated, since the energy 10 footwear can be worn on both feet and immediately

after the operation. A **removable heel** (patented system) inclined by 8° can be applied to keep the foot in talism for all rehabilitation courses where the surgeon deems it necessary. The Energy 10 footwear was subjected by the Moticon Science (independent company) of Munich, to test the **PLANTAR DISCHARGE**, comparing it with five other footwear of other brands among the most sold on the European market and was the **BEST AND MOST EFFECTIVE IN ALL THE COMPARISON TESTS.**

TOTAL OPENING THANKS
TO THE 2 ZIPS



DIFFERENT HEIGHTS OF THE
INSOLES ACCORDING TO THE NEEDS



REMOVABLE HEEL PATENTED TO KEEP
THE FOOT IN TALISM



PHASE3 PROTECTION OF FOOT HEALTH

The Energy 10 footwear has been designed for men and women who want to live their natural "vital energy". The innovative technical characteristics open the way to a new idea of well-being, health and vitality. Indeed the

line is designed to meet the needs and to prevent and manage the problems of the sensitive foot but above all, to live well.

LINE ENERGY EXTRA WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



SUITABLE FOR THE
PREVENTION OF
DIABETIC FOOT



PADDED
LINING

FOOTWEAR
→ PREVENTION
→ POST-OPERATIVE
→ REHABILITATION
→ FOR THE PROTECTION
OF FOOT HEALTH



ENERGY 10 BIO

Certification no. 2702664

Diabetic foot risk level
0 1 2 3 See page 2

Color	Natural
Sizes	35 to 43



TOTAL OPENING WITH
DOUBLE ZIP

LINE ENERGY
EXTRA WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



SUITABLE FOR THE
PREVENTION OF
DIABETIC FOOT



PADDED
LINING

ENERGY 10 BLU

Certification no. 2421679

Diabetic foot risk level

0 1 2 3

See page 2



Color

Blue

Sizes

35 to 41 (Women)
42 to 50 (Men)

ENERGY 10 NERO

Certification no. 2702666

Diabetic foot risk level

0 1 2 3

See page 2



Color

Black

Sizes

35 to 41 (Women)
42 to 50 (Men)

ENERGY PLUS

Certification no. 2702685



Color

Black

Sizes

35 to 46

Sizes Table

35-36

37-38

39-40

41-42-43

44-45-46

LINE COMFORT

Within the Comfort Line there is a collection of open footwear designed to provide maximum comfort and elegant style for both men and women. This line, which includes technical and stylish sandals and slippers, is

ideal for everyday use both indoors and outdoors, with special attention to the needs of those seeking functional footwear.

TECHNICAL FEATURES OF THE OPEN FOOTWEAR, COMFORT LINE

1. COMFORTABLE FIT

Line is designed to offer a wide and comfortable fit, able to accommodate different foot shapes, including sensitive feet or those with mild deformities.

2. ANATOMIC AND HIGH-PERFORMANCE INSOLE

The technical heart of these shoes is the anatomical insole, designed to provide optimal support:

- Arch support: Shaped to promote proper foot alignment and improve weight distribution.
- Heel cup: Stabilizes the foot and reduces impact during walking.
- Metatarsal area: Relieves pressure on the forefoot, ensuring comfort even during prolonged use.

This design makes the shoes ideal for preventing fatigue and improving the biomechanics of walking.

3. TECHNICAL AND SELF-MOLDING MATERIALS FOR THE UPPER

The upper is made of high-quality, self-molding, and breathable materials that adapt to the shape of the foot, ensuring personalized comfort. This feature is especially beneficial for those with deformities such as hallux valgus or for sensitive feet.

4. ELASTIC INSERTS ON THE UPPER

To provide additional comfort, elastic inserts are placed in the most sensitive areas (such as the big toe region) to reduce pressure and improve fitting, preventing irritation and discomfort.

5. DESIGNED FOR CUSTOM ORTHOTICS

All footwear of this line is designed to accommodate custom orthotics, providing a versatile solution that adapts to the specific needs of each customer.

6. SLIP-RESISTANT OUTSOLE

The outsole is made of high-traction materials, providing excellent safety on all surfaces. Additionally, it is lightweight but also durable, ideal for daily use, ensuring both stability and long-lasting performance.

The open footwear of the COMFORT line represents an ideal combination of technical performance, elegance, and comfort. Whether it's technical sandals for outdoor use or stylish slippers for home relaxation, this line offers a versatile and highperforming solution for those seeking COMFORT, customization, and style—without compromising foot well-being.

COMFORT LINE
WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



PADDED
LINING

TECNICA ADELE ROSSO
Certification no. 2421181

Diabetic foot risk level

0 ① ② ③

See page 2

Color Red

Sizes 35 to 43



ADJUSTABLE DOUBLE
VELCRO OPENING

TECNICA ADELE FIORI
Certification no. 2421183

Diabetic foot risk level

0 ① ② ③

See page 2

Color Flowers

Sizes 35 to 43



TECNICA ANNA FIORI
Certification no. 2421184

Diabetic foot risk level

0 ① ② ③

See page 2

Color Flowers

Sizes 35 to 43



SIDE
INSERT
STRETCH



STRETCH SIDE INSERT FOR
DIFFERENT PATHOLOGIES

COMFORT LINE
WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



PADDED
LINING

TECNICA KIARA BLU
Certification no. 2710515

Diabetic foot risk level
0 ① ② ③ See page 2

Color	Blue
Sizes	35 to 43



TECNICA ANGELA NERO
Certification no. 2710525

Diabetic foot risk level
0 ① ② ③ See page 2

Color	Black
Sizes	35 to 42



SIDE
INSERT
STRETCH



STRETCH SIDE
INSERT

TECNICA NINA
Certification no. 2710476

Diabetic foot risk level
0 ① ② ③ See page 2

Color	Black / Silver
Sizes	35 to 43



SIDE
INSERT
STRETCH



PADDED
LINING



STRETCH SIDE
INSERT

COMFORT LINE
WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



PADDED
LINING

TECNICA ALMA BLU
Certification no. 2710564

Diabetic foot risk level

0 ① ② ③

See page 2

Color Blue

Sizes 35 to 42



REAL WOOL
INSIDE

TECNICA DEVIS
Certification no. 2421194

Diabetic foot risk level

0 ① ② ③

See page 2

Color Blue

Sizes 35 to 50



ADJUSTABLE VELCRO
CLOSURE

TECNICA DYLAN
Certification no. 2421195

Diabetic foot risk level

0 ① ② ③

See page 2

Color Brown

Sizes 35 to 50



REAL WOOL INSIDE
AND ADJUSTABLE
VELCRO CLOSURE

LINE REHABILITATION

COMFORT, PROTECTION, AND VERSATILITY FOR EVERYDAY FOOT WELL-BEING

The Rehabilitation Line is designed to offer lightweight, comfortable, and technically advanced footwear, specifically intended for home use and for the needs of patients in the post-operative phase, elderly individuals, or people with chronic foot conditions. Characterized by an extra-wide fit (fitting 14) and an adjustable structure with Velcro closures, this line is ideal for managing diabetic foot. However, it also stands out for its multifunctional versatility, making it effective in the daily care of many other clinical conditions.

Indeed, the Rehabilitation Line represents a valid option for patients suffering from: hallux valgus, hammer toes, metatarsalgia, rheumatic foot, edema, post-traumatic or post-surgical deformities, and painful conditions affecting the plantar area and/or the forefoot. Each model is designed to ensure maximum facility of wear, even in the presence of bandages or swelling, and can accommodate custom orthotics up to 9 mm—providing a functional base for personalized treatments.

TECHNICAL FOCUS: OUTSOLE MODULARITY BASED ON THE CLINICAL RISK OF DIABETIC FOOT

A distinctive feature of this line is the ability to customize the outsole according to the patient's clinical profile and the rehabilitation phase:

- **FLEXIBLE OUTSOLE**, for safe and comfortable use in the home environment and in the absence of ulcers or structural issues.
- **SEMI-RIGID OUTSOLE**, designed to improve biomechanics and support the foot in moderate-risk conditions, by using an extra semi-rigid insole.
- **RIGID OUTSOLE**, for the treatment of high-risk diabetic foot and for advanced offloading and containment needs, by using an extra rigid insole.

This modularity makes the footwear highly adaptable, allowing healthcare professionals to customize the therapeutic approach based on the pathology and the environment of use.

Although designed for indoor use, the Rehabilitation Line offers a complete technical structure, featuring supportive heel counters, self-molding fabrics, and a standard breathable, antibacterial S-Air insole.

It is, in every respect, a high-performance medical shoe designed to ensure lightweight comfort, breathability, and facility of use throughout the day.

The Rehabilitation Line is a modern, customizable, and versatile therapeutic solution, capable of providing maximum support in the home setting to patients with diverse clinical needs—without compromising comfort, functionality, or protection.



Semelle intérieure Semi-Rigisupplémentaire

Certification no. 2798981

29\$



Semelle intérieure Rigisupplémentaire

Certification no. 1254508

29\$

REHABILITATION LINE EXTRA WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



SUITABLE FOR THE
PREVENTION OF
DIABETIC FOOT



ELASTIC
UPPER



PADDLED
LINING

TECNICA AGRI

Certification no. 2470730

Diabetic foot risk level

0 1 2 ③

See page 2



Color Night Blue

Sizes 35 to 50



BREATHABLE UPPER

TECNICA ANDRIA

Certification no. 2470731

Diabetic foot risk level

0 1 2 ③

See page 2



Color Green

Sizes 35 to 50



VELCRO OPENING

TECNICA FLEX

Certification no. 2702692

Diabetic foot risk level

0 1 2 ③

See page 2



Color Grey / Black

Sizes 35 to 50



REHABILITATION LINE
EXTRA WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



SUITABLE FOR THE
PREVENTION OF
DIABETIC FOOT



ELASTIC
UPPER



PADDED
LINING

TECNICA DYNAMIC AIR
Certification no. 2702570

Diabetic foot risk level
0 1 2 3 See page 2

Color	Black
Sizes	35 to 50



TECNICA SOFT
Certification no. 2702576

Diabetic foot risk level
0 1 2 3 See page 2

Color	Black
Sizes	35 to 50



STRETCH



TECNICA 3S
Certification no. 2702692

Diabetic foot risk level
0 1 2 3 See page 2

Color	Black
Sizes	35 to 50



STRETCH



REHABILITATION LINE
EXTRA WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



SUITABLE FOR THE
PREVENTION OF
DIABETIC FOOT



ELASTIC
UPPER



PADDED
LINING

TECNICA 1

Certification no. 2702539

Diabetic foot risk level

0 1 2 ③

See page 2

Color Grey

Sizes 35 to 50



TECNICA 3 AGNELLO

Certification no. 2702555

Diabetic foot risk level

0 1 2 ③

See page 2

Color Blue

Sizes 35 to 50



REAL WOOL
INSIDE

TECNICA T2 AIR

Certification no. 2702569

Diabetic foot risk level

0 1 ② ③

See page 2

Color Grey

Sizes 35 to 50



WIDE ADJUSTABLE DOUBLE
VELCRO OPENING

BAMBOO REHABILITATION LINE

MADE OF AN UPPER COVERED IN BAMBOO, A GIFT OF NATURE FOR FOOT WELL-BEING. BAMBOO IS NATURAL AND GENTLE ON THE SKIN, PREVENTS IRRITATIONS AND BAD ODORS, AND IS SUSTAINABLE AND ENVIRONMENTALLY FRIENDLY!



NATURAL ANTIBACTERIAL

Bamboo contains a natural antimicrobial agent called bamboo kun, which inhibits the growth of bacteria and fungi. This property makes bamboo ideal for direct skin contact applications, helping to prevent irritations and bad odors.



MOISTURE ABSORPTION

Bamboo has an extraordinary ability to absorb moisture, up to 3-4 times more than cotton. This helps to keep the skin dry, reducing the risk of rubbing or irritation.



BREATHABILITY AND THERMOREGULATION

Bamboo fibers are highly breathable and allow a natural temperature regulation. They keep the feet cool in summer and warm in winter, ensuring optimal comfort throughout the year.



NATURAL SOFTNESS

Bamboo fibers are incredibly soft, similar to silk, but much more durable.



HYPOALLERGENIC

Thanks to its natural composition and the absence of harsh chemicals, bamboo is an excellent choice for those with allergies or sensitive skin.



BAMBOO IS ONE OF THE MOST ENVIRONMENTALLY FRIENDLY NATURAL RESOURCES IN THE WORLD: IT GROWS RAPIDLY WITHOUT THE NEED FOR PESTICIDES OR FERTILIZERS. IT REGENERATES NATURALLY AFTER CUTTING. IT ABSORBS LARGE AMOUNTS OF CARBON DIOXIDE AND PRODUCES OXYGEN.

REHABILITATION LINE
EXTRA WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



SUITABLE FOR THE
PREVENTION OF
DIABETIC FOOT



ELASTIC
UPPER



PADDED
LINING

TECNICA 3 JAK GRIGIO
Certification no. 2755197

Diabetic foot risk level

0 1 2 ③

See page 2



Color Grey

Sizes 35 to 50



TECNICA 11 JAK GRIGIO
Certification no. 2755200

Diabetic foot risk level

0 1 2 ③

See page 2



Color Grey

Sizes 35 to 50



TECNICA 4 JAK GRIGIO
Certification no. 2755206

Diabetic foot risk level

0 1 2 ③

See page 2



Color Grey

Sizes 35 to 50



REHABILITATION LINE EXTRA EXTRA WIDE



REMOVABLE
INSOLE



EXTRA FLEXIBLE AND
SHOCK ABSORBING
OUTSOLE



PREDISPOSITION
FOR INSOLE



NO INNER
STITCHING



MEDICAL DEVICE



SUITABLE FOR THE
PREVENTION OF
DIABETIC FOOT



ELASTIC
UPPER



PADDED
LINING

XXL REHABILITATION FOOTWEAR

THIS NEW LINE OF XXL REHABILITATION FOOTWEAR WITH EXTRA-LARGE VOLUME IS DESIGNED TO PROVIDE MAXIMUM COMFORT AND SUPPORT FOR THOSE WHO NEED SHOES WITH GREAT SPACE.

Ideal for the treatment of foot pathologies, deformities, rheumatic foot, metatarsalgia, diabetic foot (risk levels 0 and 1), and post-operative rehabilitation.

MAIN FEATURES

XXL VOLUME

Designed to provide extra space and comfort, even in the presence of deformities or swelling.

ADJUSTABLE VELCRO CLOSURE

Allows a perfect and customizable fitting.

DESIGNED TO ACCOMMODATE 8 MM INSOLES

Enables the insertion of custom orthotic insoles.

COTTON LINING TREATED WITH SILVER IONS

Provides effective antibacterial treatment, helping to keep the foot dry and fresh.

SLIP-RESISTANT OUTSOLE

Ensures stability and safety during walking.

This footwear is the ideal solution for those seeking a perfect balance of comfort, protection, and functionality during the rehabilitation process or when dealing with specific foot needs.

TECNICA 13 XXL NERO

Certification no. 2755208

Diabetic foot risk level

0 1 2 3

See page 2

Color Black

Sizes 36 to 50



TECNICA 11 XXL NERO

Certification no. 2755207

Diabetic foot risk level

0 1 2 3

See page 2

Color Black

Sizes 36 to 50

